



SHOCK-RESISTANT



BREAK-PROOF



WASHABLE

PLS 001

SHOWER AND SANITIZATION PANEL



ENCLARA
MEDICAL

The PLS 001 shower panel is the perfect solution for washing patients and sanitizing equipment. It is the ideal complement to Enclara Medical shower chairs, shower trolleys and shower tubs.

BODY

The panel body is made of galvanized steel, painted with anti-corrosion epoxy powders.

SHOWER HEADS

PLS 001 is equipped with two separate handheld showerheads. Function is easily identified by different hose colour: white for showering the patient, yellow for equipment sanitization.

The handle design ensures a comfortable and ergonomic grip. The thumb lever on the handle allows to control the flow of water quickly and easily.

THERMOSTATIC MIXER

The thermostatic mixer keeps a constant water temperature throughout the shower operations. The safety button, set at 38°C, prevents scalding.

SANITIZATION

Water and disinfectant are automatically mixed. A setting valve inside the panel allows to adjust the disinfectant concentration according to the manufacturer's specifications. The sanitizer bottle housing includes a key lock to avoid improper use.

PLS 001 comes in standard version for wall installation. The pedestal is an accessory available on request.



Shower heads



Thermostatic mixer



Wall mount



Sanitizer bottle housing