



SHOCK-RESISTANT



BREAK-PROOF



WASHABLE



ANTIBACTERIAL



UV RESISTANT

CRM 001 22C

NON-MAGNETIC TROLLEY 2 + 2 DRAWERS



ENCLARA
MEDICAL

CRM 001 22C

Our CRM 001 non-magnetic trolleys are designed for the transport and distribution of medicines, dressings, instruments and various materials in all areas of the MR environment. CRM 001 trolleys are strictly made of non-metallic, non-magnetic, non-conductive materials.

For this reason, they are fully compliant with the requirements set by the EN 62570:2016 harmonized standard practice for MR SAFE marking



FRAME

The frame is entirely made of the best furniture grade PVC. The addition of titanium dioxide in the pvc compound results in a better appearance by giving it a UV-stable, glossy-white finish; it also enhances mechanical resistance and provides specific antibacterial properties.

Components are chemically welded together using a PVC-base glue, so that the frame behaves like a single part, solid, rigid and sturdy. CRM 001 is not subject to rust, scratches or dents.

DRAWERS

The drawer units, entirely in plastic, are equipped with transparent drawers with ergonomic opening and a stop system to prevent the drawers from falling. Depending on the size, each drawer has a load capacity of 15 to 20 kg.

SHELVES

The shelves are made of radiolucent HPL and are removable for easy sanitization. Each shelf has a capacity of 25 kg.

CASTORS

The trolleys are equipped with 4 large diameter swivel castors in synthetic material, 2 with total lock brake.

CRM 001 2R - CRM 001 3R - CRM 001 5C

NOTE

MR Conditional trolleys, with swivel castors in synthetic material and brass, are also available upon request, for MRI scanners up to 3 Tesla.



CRM 001 2R



CRM 001 3R



CRM 001 5C

ENCLARA MEDICAL S.r.l.

Via Liberazione, 71 - 20068 Peschiera Borromeo (Mi) - Italy - Vat Id: 10373980969
T. +39 02 94750334 - M. info@enclaramedical.it - www.enclaramedical.it