



SHOCK-RESISTANT



BREAK-PROOF



WASHABLE



ANTIBACTERIAL



UV RESISTANT

# CRM 001 3R

NON-MAGNETIC TROLLEY 3 SHELVES



**ENCLARA**  
**MEDICAL**

# CRM 001 3R

Our CRM 001 non-magnetic trolleys are designed for the transport and distribution of medicines, dressings, instruments and various materials in all areas of the MR environment. CRM 001 trolleys are strictly made of non-metallic, non-magnetic, non-conductive materials.

**For this reason, they are fully compliant with the requirements set by the EN 62570:2016 harmonized standard practice for MR SAFE marking**



## FRAME

The frame is entirely made of the best furniture grade PVC. The addition of titanium dioxide in the pvc compound results in a better appearance by giving it a UV-stable, glossy-white finish; it also enhances mechanical resistance and provides specific antibacterial properties.

Components are chemically welded together using a PVC-base glue, so that the frame behaves like a single part, solid, rigid and sturdy. CRM 001 is not subject to rust, scratches or dents.



## SHELVES

The shelves are made of radiolucent HPL and are removable for easy sanitization. Each shelf has a capacity of 25 kg.

## CASTORS

The trolleys are equipped with 4 large diameter swivel castors in synthetic material, 2 with total lock brake.

CRM 001 2R - CRM 001 5C - CRM 001 22C

## NOTE

MR Conditional trolleys, with swivel castors in synthetic material and brass, are also available upon request, for MRI scanners up to 3 Tesla.



CRM 001 2R



CRM 001 5C



CRM 001 22C

ENCLARA MEDICAL S.r.l.

Via Liberazione, 71 - 20068 Peschiera Borromeo (Mi) - Italy - P.iva: 10373980969  
T. +39 02 94750334 - M. info@enclaramedical.it - www.enclaramedical.it